


Product Name: Recombinant Human SEMA5A Protein
Catalog #: PEH6219



Summary

Name	SEMA5A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q13591
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TSWSQCSRDC SRGIRNRKRV CNNPEPKYGG MPCLGPSLEY QECNILPCPV DGVWSCWSPW TKCSATCGGG HYMRTRSCSN PAPAYGGDIC LGLHTEEALC NTQPCPESWS EWSDWSECEA SGVQVRARQC ILLFPMGSQC SGNTTESRPC VFDSNFIPEV SVARSSSVEE

image



Background

Alternative Names Semaphorin-5A (Semaphorin-F) (Sema F)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SEMA5A Protein
Catalog #: PEH6219

