


Product Name: Recombinant Human SPIN2B Protein
Catalog #: PEH6220



Summary

Name	SPI2B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BPZ2
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MKTPNAQEAE GQQTAAAGR ATGSANMTKK KVSQKKQRGR PSSQPRRNIV GCRISHGWKE GDEPITQWKG TVLDQVPINP SLYLVKYDGI DCVYGLELHR DERVLSLKIL SDRVASSHIS DANLANTIIG KAVEHMFEGE HGSKDEWRGM VLAQAPIMKA WFYITYEKDP VLYMYQLLDD YKEGDLRIMP ESSESPTER EPGGVVDGLI GKHVEYTKED GSKRIGMVIH QVEAKPSVYF IKFDDDFHIY VYDLVKKS

image



Background

Alternative Names	Spindlin-2B (Spindlin-like protein 2B) (SPIN-2) (SPIN-2B)
--------------------------	---

Product Name: Recombinant Human SPIN2B Protein
Catalog #: PEH6220



Note

For Research Use Only , Not for Diagnostic Use