


Product Name: Recombinant Human LY6G6D Protein
Catalog #: PEH6231



Summary

Name	LY66D_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95868
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NRMRCYNCGG SPSSSCKEAV TTCGEGRPQP GLEQIKLPGN PPVTLIHQHP ACVAAHHCNQ VETESVGDVT YPAHRDCYLG DLCNS

image



Background

Alternative Names	Lymphocyte antigen 6 complex locus protein G6d (Protein Ly6-D) (Megakaryocyte-enhanced gene transcript 1 protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use