


Product Name: Recombinant Human LRFN3 Protein
Catalog #: PEH6237



Summary

Name	LRFN3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9BTN0
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LHCNCELVWL RRLAREDDLE ACASPPALGG RYFWAVGEEE FVCEPPVVTH RSPPLAVPAG RPAALRCRAV GDPEPRVRWV SPQGRLGNS SRARAFPNGT LELLVTEPGD GGIFTCIAAN AAGEATAAVE LTVGPPPPPQ LANSTSCDPP RDGDPDALTP

image



Background

Alternative Names Leucine-rich repeat and fibronectin type-III domain-containing protein 3 (Synaptic adhesion-like molecule 4)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human LRFN3 Protein
Catalog #: PEH6237

