


Product Name: Recombinant Human NDUB6 Protein
Catalog #: PEH6239



Summary

Name	NDUB6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95139
Host	E.coli
Species	Human
Molecular Mass	9 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RLQQLRELRR RWLKDQELSP REPVLPQKM GPMEKFWNKF LENKSPWRKM

image



Background

Alternative Names NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 6 (Complex I-B17) (CI-B17) (NADH-ubiquinone oxidoreductase B17 subunit)

Note

For Research Use Only , Not for Diagnostic Use