


Product Name: Recombinant Human TBX1 Protein
Catalog #: PEH6243



Summary

Name	TBX1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43435
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PFAPAAGAAT SAAAEPEGPG ASCAAAKAP VKKNAKVAGV SVQLEMKALW DEFNQLGTEM IVTKAGRRMF PTFQVKLFGM DPMADYMLLM DFVPVDDKRY RYAFHSSSWL VAGKADPATP GRVHYHPDSP AKGAQWMKQI VSFDKLKLTN NLLDDNGHII LNSMHRYQPR FHVVYVDPRK DSEKYAEENF KTFVFEETRF TAVTAYQNHR ITQLKIASNP FAKGFRDCDP EDWPRNHRPG ALPLMSAFAR SRNPVASPTQ PSGTEKGGHV LKDKEVKAET SRNTPEREVE LLRDAGGCVN LGLPCPAECQ PFNTQGLVAG

image



Background

Alternative Names	T-box transcription factor TBX1 (T-box protein 1) (Testis-specific T-box protein)
--------------------------	---

Product Name: Recombinant Human TBX1 Protein
Catalog #: PEH6243



Note

For Research Use Only , Not for Diagnostic Use