


Product Name: Recombinant Human TBX4 Protein
Catalog #: PEH6244



Summary

Name	TBX4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P57082
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEAFRAPGPA LGEASAANAP EPALAAPGLS GAALGSPPGP GADVVA EQTNIKVG LHEKELWKKF HEAGTEMIIT KAGRRMFPSY KVKVTGMNPK TKYILLIDIV PADDHRYKFC DNKWMVAGKA EPAMPGRLYV HPDSPATGAH WMRQLVSFQK LKLTNNHLDP FGHIILNSMH KYQPRHLIVK ADENNAFGSK NTAFCTHVFP ETSFISVTSY QNHKITQLKI ENNPFKGFGR GSDDSDLRVA

image



Background

Alternative Names T-box transcription factor TBX4 (T-box protein 4)

Note

Product Name: Recombinant Human TBX4 Protein
Catalog #: PEH6244



For Research Use Only , Not for Diagnostic Use