


Product Name: Recombinant Human TENM4 Protein
Catalog #: PEH6248



Summary

Name	TEN4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6N022
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EITEDTASSW PVPTDVS LYP SGGTGLET PD RKGKGTTEGK PSSFFPEDSF IDSGEIDVGR RASQKIPPGT FWRSQVFIDH PVHLKFNVSL GKAALVGIYG RKGLPPSHTQ FDFVELLDGR RLLTQEARSL EGTPRQSRGT VPPSSHETGF IQYLDSGIWH LAFYNDGKES

image



Background

Alternative Names	Teneurin-4 (Ten-4) (Protein Odd Oz/ten-m homolog 4) (Tenascin-M4) (Ten-m4) (Teneurin transmembrane protein 4)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TENM4 Protein
Catalog #: PEH6248

