



Product Name: Recombinant Human ZNF431 Protein
Catalog #: PEH6264



Summary

Name	ZN431_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TF32
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SYKCEECGKG FNWSSTLTKH KRIHTGEKPY KCEVCGKAFN ESSNLTHKM IHTGEKPYKC EECGKAFNRS PQLTAHKIIH TGEKPYKCEE CGKAFSQSSI LTTHKRIHTG EKPYKCEECG KAFNRSSNLT KHKIIHTGEK SYKCEECGKA FNQSSTLTKH

image



Background

Alternative Names Zinc finger protein 431

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZNF431 Protein
Catalog #: PEH6264

