


**Product Name: Recombinant Human ADCYAP1 Protein**  
**Catalog #: PEH6270**



## Summary

|                              |   |
|------------------------------|---|
| <b>Name</b>                  | PACA_HUMAN  |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE   |
| <b>Endotoxin level</b>       | Please contact with the lab for this information  |
| <b>Tag</b>                   | C-His   |
| <b>Accession #</b>           | P18509  |
| <b>Host</b>                  | E.coli  |
| <b>Species</b>               | Human   |
| <b>Molecular Mass</b>        | 11 kDa  |
| <b>Buffer</b>                | Lyophilized from Phosphate buffered saline  |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.  |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| <b>Sequence</b>              | RDVAHGILNE AYRKVLDQLS AGKHLQSLVA RGVGGSLGGG AGDDAEPLSK<br>RHSDGIFTDS YSRYRKQMAV KKYLA AVL GK RYKQ RVKNKG  |

## image



## Background

|                          |   |
|--------------------------|---|
| <b>Alternative Names</b> | Pituitary adenylate cyclase-activating polypeptide (PACAP) [Cleaved into: PACAP-related peptide (PRP-48); Pituitary adenylate cyclase-activating polypeptide 27 (PACAP-27) (PACAP27); Pituitary adenylate cyclase-activating polypeptide 38 (PACAP-38) (PACAP38)] |
|--------------------------|---|

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human ADCYAP1 Protein**  
**Catalog #: PEH6270**

