


Product Name: Recombinant Human UHMK1 Protein
Catalog #: PEH6279



Summary

Name	UHMK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TAS1
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DDYLENEEY EDVVEDVKEE CQKYGPVSL LVPKENPGRG QVFVEYANAG DSKAAQKLLT GRMFDGKFVV

image



Background

Alternative Names Serine/threonine-protein kinase Kist (EC 2.7.11.1) (Kinase interacting with stathmin) (PAM COOH-terminal interactor protein 2) (P-CIP2) (U2AF homology motif kinase 1)

Note

For Research Use Only , Not for Diagnostic Use