


Product Name: Recombinant Human UBR1 Protein
Catalog #: PEH6287



Summary

Name	UBR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8I WV7
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	THDKEKAERK RKAEARLHR QKIM AQMSAL QKNFIETHKL MYDNTSEMPG KEDSIM EES TPAVSDYSRI ALGPKRGPSV TEKEVLTCL CQEEQEVKIE

image



Background

Alternative Names E3 ubiquitin-protein ligase UBR1 (EC 2.3.2.27) (N-recognin-1) (RING-type E3 ubiquitin transferase UBR1) (Ubiquitin-protein ligase E3-alpha-1) (Ubiquitin-protein ligase E3-alpha-I)

Note

For Research Use Only , Not for Diagnostic Use