



Product Name: Recombinant Human VWC2L Protein
Catalog #: PEH6291



Summary

| | |
|------------------------------|---|
| Name | VWC2L_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | N-SUMO & C-His |
| Accession # | B2RUY7 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 25 kDa |
| Buffer | Lyophilized from Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | KNYKILEEFK PSPCEWCRCE PSNEVHCVVA DCAVPECVNP VYEPEQCCPV CKNGPNCFAG TTIPAGIEV KVDECNICHC HNGDWWKPAQ CSKRECQGKQ |

image



Background

Alternative Names von Willebrand factor C domain-containing protein 2-like (Brorin-like)

Note

For Research Use Only , Not for Diagnostic Use