


Product Name: Recombinant Human TPSG1 Protein
Catalog #: PEH6297



Summary

Name	TRYG1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9NRR2
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VSDAGGRIVG GHAAPAGAWP WQASRLRRV HVCGGSLSP QWVLTAHCF SGSLNSSDYQ VHLGELEITL SPHFSTVRQI ILHSSPSGQP GTSGDIALVE LSVPVTLSSR ILPVCLPEAS DDFCPGIRCW VTGWGYTREG EPLPPYSLR EVKVSVDTE TCRDYPGPG GSILQPDMLC ARGPGDACQD DSGGPLVCQV NGAWVQAGTV SWGEGCGRPN RPGVYTRVPA YVNWIRRHIT ASGGSESGYP

image



Background

Alternative Names	Tryptase gamma (EC 3.4.21.-) (Serine protease 31) (Transmembrane tryptase) [Cleaved into: Tryptase gamma light chain; Tryptase gamma heavy chain]
--------------------------	---

Product Name: Recombinant Human TPSG1 Protein
Catalog #: PEH6297



Note

For Research Use Only , Not for Diagnostic Use