

Product Name: Recombinant Human IGSF9B Protein
Catalog #: PEH6303



Summary

Name	TUTLB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UPX0
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RAPSESSDDQ GQPAKRMLS PTREKELSLY KKTKRAISSK KYSVAKAEAE AEATTPIELI SRGPDGRFVM DPAEMEPLSK SRRIEGFPFA EETDMYPEFR

image



Background

Alternative Names Protein turtle homolog B (Immunoglobulin superfamily member 9B) (IgSF9B)

Note

For Research Use Only , Not for Diagnostic Use