


Product Name: Recombinant Human ZNF430 Protein
Catalog #: PEH6328



Summary

Name	ZN430_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H8G1
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LKEASGCPGA DRNLLVYSFY EKGPLTFRDV AIEFSLEEWQ CLDTAQQDLY RKVMLENYRN LVFLAGIAVS KPDLITCLEQ GKEPWNMKRH AMVDQPPVTY SHFAQDLWPE QGIKDSFQEV ILRRYGKCGH EDLQLRTGCK SVDECNLHKE

image



Background

Alternative Names Zinc finger protein 430

Note

For Research Use Only , Not for Diagnostic Use