


Product Name: Recombinant Human ZNF648 Protein
Catalog #: PEH6329



Summary

Name	ZN648_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5T619
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GEASPLSSLT EEAHDTQMLS MNLESDDDEDG GEAEKEGTAD PVACPRGSSP VTHENPDLPW PHPLGKEEEK FSDSSSAGGM GQKPVEMSGK ASWSRDVTKI NETQGSPGAS RALGSLPSGL AHKLLGQMGP LGDRLPAGDD GYSGANQDAV LDVPPSFPSN GKYLCAHKSV DTSAGNSSLL CFPRPGSNWD LPTQETHTPA QASATPASLA AAVLAKARNS RKVQNQAGRR EGGEAEARPY RCLRGGRAFQ

image



Background

Alternative Names Zinc finger protein 648

Note

Product Name: Recombinant Human ZNF648 Protein
Catalog #: PEH6329



For Research Use Only , Not for Diagnostic Use