


Product Name: Recombinant Human USP19 Protein
Catalog #: PEH6334



Summary

Name	UBP19_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O94966
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RRGPPGLEDT TSKKKQKDRA NQESKDGDPK KETGSRVVAQ AGLEPLASGD PSASASHAAG ITGSRHRTLL FFPSSSGSAS TPQEEQTKEG

image



Background

Alternative Names	Ubiquitin carboxyl-terminal hydrolase 19 (EC 3.4.19.12) (Deubiquitinating enzyme 19) (Ubiquitin thioesterase 19) (Ubiquitin-specific-processing protease 19) (Zinc finger MYND domain-containing protein 9)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use