


Product Name: Recombinant Human TMEM79 Protein
Catalog #: PEH6341



Summary

Name	TMM79_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BSE2
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EVKRSDSPEK SSPQALVPNG RQPEGEGGAE SPGAESLRVG SSAGSPTAIE GAEDGLDSTV SEAATLPWGT GPQPSAPFPD PPGWRDIEPE PPESEPLTKL EELPEDDANL LPEKAARAFV PIDLQCIERQ PQEDLIVRCE AGEGECRTFM PPRVTHPDPT ERKWAEAVVR PPGCSCGGCG SCGDREWLRA

image



Background

Alternative Names Transmembrane protein 79 (Matrin)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TMEM79 Protein
Catalog #: PEH6341

