


Product Name: Recombinant Human TGFBR1 Protein
Catalog #: PEH6342



Summary

Name	TGFR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P36897
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CTKDNFTCVT DGLCFVSVTE TTDKVIHNSM CIAEIDLIPR DRPFVCAPSS KTGSVTTTYC CNQDHCNKIE LPTTVKSSPG

image



Background

Alternative Names	TGF-beta receptor type-1 (TGFR-1) (EC 2.7.11.30) (Activin A receptor type II-like protein kinase of 53kD) (Activin receptor-like kinase 5) (ALK-5) (ALK5) (Serine/threonine-protein kinase receptor R4) (SKR4) (TGF-beta type I receptor) (Transforming growth factor-beta receptor type I) (TGF-beta receptor type I) (TbetaR-I)
--------------------------	---

Note

Product Name: Recombinant Human TGFBR1 Protein
Catalog #: PEH6342



For Research Use Only , Not for Diagnostic Use