


Product Name: Recombinant Human PRSS21 Protein
Catalog #: PEH6351



Summary

Name	TEST_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9Y6M0
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IVGGEDAELG RWPWQGSRL WDSHVCVSL LSHRWALTAA HCFETSDLS DPSGWMVQFG QLTSMPFSWS LQAYYTRYFV SNIYLSRYL GNSPYDIALV KLSAPVTYTK HIQPICLQAS TFEFENRTDC WVTGWGYIKE DEALPSPTL QEVQVAIINN SMCNHLFLKY SFRKDIFGDM VCAGNAQGGK DACFGDSGGP LACNKNGLWY QIGVVSWGVG CGRPNRPGVY TNISHHFEWI QKLMAQS

image



Background

Alternative Names	Testisin (EC 3.4.21.-) (Eosinophil serine protease 1) (ESP-1) (Serine protease 21)
--------------------------	--

Note

Product Name: Recombinant Human PRSS21 Protein
Catalog #: PEH6351



For Research Use Only , Not for Diagnostic Use