


Product Name: Recombinant Human VSX2 Protein
Catalog #: PEH6353



Summary

Name	VSX2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P58304
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ASGPLDTSQT ASSDSEDVSS SDRKMSKSAL NQTKKRKRR HRTIFTSYQL EELEKAFNEA HYPDVYAREM LAMKTELPED RIQVWFQNRK AKWRKREKCW GRSSVMAEYG LYGAMVRHSI PLPESILKSA KDGIMDSCAP WLLGMHKKSL EAAAESGRKP EGERQALPKL DKMEQDERGP DAQAAISQEE LRENSIAVLR AKAQEHSTKV LGTVSGPDSL ARSTEKPEEE EAMDEDRPAE

image



Background

Alternative Names Visual system homeobox 2 (Ceh-10 homeodomain-containing homolog) (Homeobox protein CHX10)

Note

Product Name: Recombinant Human VSX2 Protein
Catalog #: PEH6353



For Research Use Only , Not for Diagnostic Use