


Summary

Name	VPS16_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H269
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDVARAINQK LGDTPGVSYS DIAARAYGCG RTELAIKLE YEPRSGEQVP LLLKMKRSKL ALSKAIESGD TDLVFTVLLH LKNELNRGDF FMTLRNQPMA LSLYRQFCKH QELETLDKLY NQDDNHQELG SFHIRASYAA EERIEGRVAA LQTAADAFYK AKNEFAAKAT EDQMRLRLQ RRLEDELGGQ FLDLSLHDTV TTLILGGHNK RAEQLARDFR IPDKRLWWLK LTALADLEDW EELEKFSKSK KSPIGYLPFV EICMKQH NKY EAKKYASRVG

image



Background

Alternative Names Vacuolar protein sorting-associated protein 16 homolog (hVPS16)

Product Name: Recombinant Human VPS16 Protein
Catalog #: PEH6359



Note

For Research Use Only , Not for Diagnostic Use