


Product Name: Recombinant Human TMX1 Protein
Catalog #: PEH6362



Summary

Name	TMX1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H3N1
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VRVITDENWR ELLEGDWMIE FYAPWCPACQ NLQPEWESFA EWGEDLEVNI AKVDVTEQPG LSGRFIITAL PTIYHCKDGE FRRYQGPRTK KDFINFISDK

image



Background

Alternative Names	Thioredoxin-related transmembrane protein 1 (Thioredoxin domain-containing protein 1) (Transmembrane Trx-related protein)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use