


Product Name: Recombinant Human SLC30A9 Protein
Catalog #: PEH6364



Summary

Name	ZNT9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6PML9
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GRSIQPEQVQ RLTELENDP SVRAIHDVKA TDLGLGKVRV KAEVDFDGRV VTRSYLEKQD FDQMLQEIQE VKTPEELET F MLKHGENIID TLGAEVDRLE KELKKRNPEV

image



Background

Alternative Names Zinc transporter 9 (ZnT-9) (Human embryonic lung protein) (HuEL) (Solute carrier family 30 member 9)

Note

For Research Use Only , Not for Diagnostic Use