


Product Name: Recombinant Human PRPS2 Protein
Catalog #: PEH6366



Summary

Name	PRPS2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P11908
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SHQDLSQRVA DRLGLELGKV VTKKFSNQET SVEIGESVRG EDVYIIQSGC GEINDNLMEL LIMINACKIA SSSRVTAVIP CFPYARQDKK

image



Background

Alternative Names	Ribose-phosphate pyrophosphokinase 2 (EC 2.7.6.1) (PPRibP) (Phosphoribosyl pyrophosphate synthase II) (PRS-II)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use