



Product Name: Recombinant Human ZNF22 Protein
Catalog #: PEH6368



Summary

| | |
|------------------------------|--|
| Name | ZNF22_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | P17026 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 14 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | YKDECGESF KQSSNLIQH Q RIHTGEKPYQ CDECGRCFSQ SSHLIQHQR T HTGEKPYQCS ECGKCF SQSS HLRQHMKVHK EEKPRKTRGK NIRVKTHLPS |

image



Background

Alternative Names Zinc finger protein 22 (Zinc finger protein KOX15) (Zinc finger protein Krox-26)

Note

For Research Use Only , Not for Diagnostic Use