


Product Name: Recombinant Human NFIX Protein
Catalog #: PEH6377



Summary

Name	NFIX_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14938
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PSYYNINQVT LGRRSITSP STSTTKRPKS IDDSEMEPV DDVFYPGTGR SPAAGSSQSS GWPNDVDAGP ASLKKSGKLD FCSALSSQGS

image



Background

Alternative Names Nuclear factor 1 X-type (NF1-X) (Nuclear factor 1/X) (CCAAT-box-binding transcription factor) (CTF) (Nuclear factor I/X) (NF-I/X) (NFI-X) (TGGCA-binding protein)

Note

For Research Use Only , Not for Diagnostic Use