


Product Name: Recombinant Human SNRNP35 Protein
Catalog #: PEH6385



Summary

Name	U1SBP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16560
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VKNDLYREGK RERRERSR ERHWDSRTRD RDHDRGREKR WQEREPTRVW PDNDWERERD FRDDRIKGRE

image



Background

Alternative Names	U11/U12 small nuclear ribonucleoprotein 35 kDa protein (U11/U12 snRNP 35 kDa protein) (U11/U12-35K) (Protein HM-1) (U1 snRNP-binding protein homolog)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use