


Product Name: Recombinant Human ZNF529 Protein
Catalog #: PEH6399



Summary

Name	ZN529_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6P280
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LRDVVINFSQ EEWEYLDSAQ RNLYWDVMME NYSNLLSLDL ESRNETKHLS VGKDIIQNTG SQWEVMESSK LCGLEGSIFR NDWQSKSKID LQGPEVGYFS QMKIISENVP SYKTHESLTL PRRTHDSEKP YEYKEYEKFV SCDLEFDEYQ KIHTGGKKNYE CNQCWKTFGI

image



Background

Alternative Names Zinc finger protein 529

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZNF529 Protein
Catalog #: PEH6399

