


Product Name: Recombinant Human ZNF14 Protein
Catalog #: PEH6401



Summary

Name	ZNF14_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P17017
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VNFTLEEWAL LDSSQKKLYE DVMQETFKNL VCLGKKWEDQ DIEDDHRNQG KNRRCHMVER LCESRRGSKC

image



Background

Alternative Names Zinc finger protein 14 (Gonadotropin-inducible ovary transcription repressor 4) (GIOT-4) (Zinc finger protein KOX6)

Note

For Research Use Only , Not for Diagnostic Use