


Product Name: Recombinant Human ZNF683 Protein
Catalog #: PEH6404



Summary

Name	ZN683_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IZ20
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PGLQASSTDD KKFTVKYPQN KDKLGKQPER AGEGAPCPAF SSHNSSSPPP LQNRKSPSPL AFCPCPPVNS

image



Background

Alternative Names Tissue-resident T-cell transcription regulator protein ZNF683 (Homolog of Blimp-1 in T-cell) (Hobit) (Zinc finger protein 683)

Note

For Research Use Only , Not for Diagnostic Use