


Summary

Name	ZSC29_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IWIY8
Host	E.coli
Species	Human
Molecular Mass	66 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GTNSETFRQR FRRFHYQEVA GPREAFSQLW ELCCRWLRPE VRTKEQIVEL LVLEQFLTVL PGEIQNWVQE QCPENGEEAV TLVEDLEREP GRPRSSVTVS VKGQEVRLK MTPPKSSQEL LSVRQESVEP QPRGVPKKER ARSPDLGPQE QMNPKEKLP FQRSGLPFPK SGVVSRLQEQ EPWIPDLLGS KEKELPSGSH IGDRRVHADL LPSKKDRRSW VEQDHWSFED EKVAGVHWGY EETRLLAIL SQTEFYEARL NCHRNSQVYG AVAERLREYG FLRTLEQCRT KFKGLQKSYR KVKSGHPPET CPFFEMEAL MSAQVIALPS NGLEAAASHS GLVGSDAETE EPGQRGWQHE EGAEAEVAQE SDSDDMDLEA TPQDPNSAAP VVFRSPGGVH WGYEETKTYL AILSETQFYE ALRNCHRNSQ LYGAVAERLW EYGFLRTPEQ CRTKFKSLQT SYRKVKNQQA PETCPFFEEM DALVSVRVAA PPNDGQEETA SCPVQGTSEA EAQKQAEAD EATEEDSDDD EEDTEIPPGA VITRAPVLFQ SPRGFEAGFE NEDNSKRDIS

image



Product Name: Recombinant Human ZSCAN29 Protein
Catalog #: PEH6405



Background

Alternative Names Zinc finger and SCAN domain-containing protein 29 (Zinc finger protein 690)

Note

For Research Use Only , Not for Diagnostic Use