


Product Name: Recombinant Human ZNF296 Protein
Catalog #: PEH6406



Summary

Name	ZN296_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WUU4
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSKLNRRHKKT HRQVPPQSPL MADTSQEQAS AAPPEPAVHA AAPTSTLPCS GGEGAGAAAT AGVQEPGAPG SGAQAGPGGD TWGAIITTEQR TDPANSQKAS PKKMPKSGGK SRGPGGSCEF CGKHFTNSSN

image



Background

Alternative Names Zinc finger protein 296 (ZFP296) (Zinc finger protein 342)

Note

For Research Use Only , Not for Diagnostic Use