


Product Name: Recombinant Human ZSCAN10 Protein
Catalog #: PEH6409



Summary

Name	ZSC10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96SZ4
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DVTLEEKGCA SQVPSHSPKK ELPAEEPSVL GPSDEPPRPQ PRAAQPAEPG QWRLPPSSKQ PLSPGPQKTF QALQEESPQG PSPWPPEESSR DQELAAVLEC LTFEDVPENK AWPAPHLGFG SRTPDKKEEFK QEEPKGAAWP TPILAESQAD SPGVPGEPKA QSLGRGAAAS GPGEDGSLLG SSEILEVKVA EGVPEPNPEL QFICADCGVS FPQLSRLKAH QLRSHAPAGRS FLCLCCGKSF GRSSILKLHM

image



Background

Alternative Names Zinc finger and SCAN domain-containing protein 10 (Zinc finger protein 206)

Note

Product Name: Recombinant Human ZSCAN10 Protein
Catalog #: PEH6409



For Research Use Only , Not for Diagnostic Use