


Product Name: Recombinant Human ZSWIM7 Protein
Catalog #: PEH6411



Summary

Name	ZSWIM7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q19AV6
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAVVLPAVVE ELLSEMAAAV QESARIPDEY LLSLKFLFGS SATQALDLVD RQSITLISSP SGRRVYQVLG SSSKTYTCLA SCHYCSCPAF AFSVLRKSDS ILCKHLLAVY LSQVMRTCQQ LSVSDKQLTD ILLMEKKQEA

image



Background

Alternative Names	Zinc finger SWIM domain-containing protein 7 (SWIM domain-containing and Srs2-interacting protein 1 homolog) (SWIM-type zinc finger domain-containing protein 7)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZSWIM7 Protein
Catalog #: PEH6411

