


Product Name: Recombinant Human LIN7C Protein
Catalog #: PEH6417



Summary

Name	LIN7C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NUP9
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PRVVELPKTE EGLGFNIMGG KEQNSPIYIS RIIPGGIADR HGGLKRGDQL LSVNGVSVEG EHHEKAVELL KAAQGKVKLV VRYTPKVL EE MESRF EKMRS

image



Background

Alternative Names Protein lin-7 homolog C (Lin-7C) (Mammalian lin-seven protein 3) (MALS-3) (Vertebrate lin-7 homolog 3) (Veli-3)

Note

For Research Use Only , Not for Diagnostic Use