


Product Name: Recombinant Human TRPM2 Protein
Catalog #: PEH6418



Summary

Name	TRPM2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O94759
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NYLQNRQFQQ KQRPEQKIED ISNKVDAMVD LLDLDPLKRS GSMEQRLASL EEQVAQTAQA LHWIVRTLRA SGFSSEADVP TLASQKAAEE PDAEPGGRKK

image



Background

Alternative Names	Transient receptor potential cation channel subfamily M member 2 (Estrogen-responsive element-associated gene 1 protein) (Long transient receptor potential channel 2) (LTrpC-2) (LTrpC2) (Transient receptor potential channel 7) (TrpC7) (Transient receptor potential melastatin 2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TRPM2 Protein
Catalog #: PEH6418

