


Product Name: Recombinant Human OLFM3 Protein
Catalog #: PEH6419



Summary

Name	NOE3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96PB7
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DTLTQISPKE GWQVYSSAQD PDGRCICTVV APEQNLCSR D AKSRQLRQLL EKVQNMSQSI EVLNLRTQRD FQYVLKMETQ MKGLKAKFRQ IEDDRKTLMT KHFQELKEKM DELLPLIPVL EQYKTDAKLI TQFKEIRNL SAVLTGIQEE IGAYDYEELH QRVLSLETRL

image



Background

Alternative Names Noelin-3 (Olfactomedin-3) (Optimedin)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human OLFM3 Protein
Catalog #: PEH6419

