


Product Name: Recombinant Human PIGS Protein
Catalog #: PEH6428



Summary

Name	PIGS_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96S52
Host	E.coli
Species	Human
Molecular Mass	47 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ALDHEEEALS SGSVQEAEM LDEPQEQAEG SLTVYVISEH SLLPQDMMS YIGPKRTAVV RGIMHREAFN IIGRRIVQVA QAMSLTEDVL AAALADHLPE DKWSAEKRRP LKSSLGYEIT FLLNPDPS HDVYWDIEGA VRRYVQPFLN ALGAAGNFSV DSQILYYAML GVNPRFDSAS SSYYLDMHSL PHVINPVESR LGSSAASLYP VLNFLLYVPE LAHSPLYIQD KDGAPVATNA FHSPRWGGIM VYNVDSKTYN ASVLPVRVEV DMVRVMEVFL AQLRLLFGIA QPQLPPKCLL SGPTSEGLMT WELDRLLWAR SVENLATATT TLTSLAQLLG KISNIVIKDD VASEVYKAVA AVQKSAEELA SGHLASAFVA SQEAVTSSEL AFFDPSLLHL

image



Background

Alternative Names GPI transamidase component PIG-S (Phosphatidylinositol-glycan

Product Name: Recombinant Human PIGS Protein
Catalog #: PEH6428



biosynthesis class S protein)

Note

For Research Use Only , Not for Diagnostic Use