



**Product Name: Recombinant Human PLAAT4 Protein**  
**Catalog #: PEH6430**



## Summary

<b>Name</b>	PLAT4_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9UL19
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	16 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	GDLIEIFRLG YEHWALYIGD GYVIHLAPPS EYPGAGSSSV FSVLSNSAEV KRERLEDVVG GCCYRVNNSL DHEYQPRPVE VISSAKEMV GQKMKYSIVS RNCEHFVTQL RYGKSRCKQV

## image



## Background

<b>Alternative Names</b>	Phospholipase A and acyltransferase 4 (EC 2.3.1.-) (EC 3.1.1.32) (EC 3.1.1.4) (HRAS-like suppressor 4) (HRSL4) (RAR-responsive protein TIG3) (Retinoic acid receptor responder protein 3) (Retinoid-inducible gene 1 protein) (Tazarotene-induced gene 3 protein)
--------------------------	---

## Note

**Product Name: Recombinant Human PLAAT4 Protein**  
**Catalog #: PEH6430**



---

For Research Use Only , Not for Diagnostic Use