



Product Name: Recombinant Human TECTB Protein
Catalog #: PEH6432



Summary

Name	TECTB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q96PL2
Host	E.coli
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	APNKADVILV FCYPKTIITK IPECPYGWEV HQLALGGLCY NGVHEGGYYQ FVIPDLSPKN KSYCGTQSEY KPPIYHFYSH IVSNDTTVIV KNQPVNYSFS CTYHSTYLVN QAAFDQRVAT VHVKNQSMGT FESQLSLNFY TNAKFSIKKE APFVLEAEI GSDLFAGVEA KGLSIRFKVV LNSCWATPSA DFMYPQLWQL INKGCPTDET VLVHENGGRDH RATFQFNAFR FQNIPKLSKV WLHCETFICD SEKLSCPVTC

image



Background

Alternative Names Beta-tectorin

Product Name: Recombinant Human TECTB Protein
Catalog #: PEH6432



Note

For Research Use Only , Not for Diagnostic Use