


Product Name: Recombinant Human ZNF234 Protein
Catalog #: PEH6474



Summary

Name	ZN234_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14588
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KDVAVVFTEE ELGLDPVQR NLYQDVMLEN FRNLLSVGHH PFKHDFVLE KEKKLDIMKT ATQRKGKSAD KIQSEVETVP EAGRHEELYW GQIWKQIASD LIKYEDSMIS ISRFPRQGD L SCQVRAGLYT THTGQKFYQC DEYKKSFTDV FNFDLHQQ LH SGEKSHTCDE

image



Background

Alternative Names Zinc finger protein 234 (Zinc finger protein 269) (Zinc finger protein HZF4)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZNF234 Protein
Catalog #: PEH6474

