


Product Name: Recombinant Human ZNF248 Protein
Catalog #: PEH6475



Summary

Name	ZN248_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NDW4
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KDVCVDFTQE EWYLLDPAQK ILYRDVILEN YSNLVSVMGYC ITKPEVIFKI EQGEPPWILE KGFPSQCHPE RKWKVDDVLE SSQENEDDHF WELLFHNNKT VSVENGDRGS KTFNLGTDVP SLRNYPYKIC DSCMNLLKNI SGLIISKKNC SRKKPDEFNV CEKLLDIRH EKIPIGEKSY KYDQKRNAIN YHQDLSQPSF GQSFEYSKNG QGFHDEAAFF TNKRSQIGET VCKYNECGRT FIESLKLNIS QRPHLEMEPY

image



Background

Alternative Names	Zinc finger protein 248
--------------------------	-------------------------

Product Name: Recombinant Human ZNF248 Protein
Catalog #: PEH6475



Note

For Research Use Only , Not for Diagnostic Use