


Product Name: Recombinant Human ZNF749 Protein
Catalog #: PEH6476



Summary

Name	ZN749_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43361
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDVAIYFSQE EWGILNDAQR HLHSNVMLEN FALLSSVGCW HGAKDEEVPS KQCVSVRVLQ VTIPKPALST LKAQPCKMCS SILKDILHLA EHDGTHPEQG LYTCAA EHDL HQKEQIREKL TRSDEWRPSF VNHS AHVGER NFTCTQGGKD FTASSDLLQQ QVLNSGWKLY RDTQDGEAFQ GEQNDFNSSQ G GKDFCHQHG LFEHQKTHNG

image



Background

Alternative Names Zinc finger protein 749

Note

Product Name: Recombinant Human ZNF749 Protein
Catalog #: PEH6476



For Research Use Only , Not for Diagnostic Use