


Product Name: Recombinant Human ZNF275 Protein
Catalog #: PEH6484



Summary

Name	ZN275_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NSD4
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AKYSESTSAT RHQMKGEDAQ PQEMASTSFP RASGPSPEFR QHGDSGKRG SPQNLPIEHH FACKECGDTF RLKVLVQHQ RVHSEKQWE CGDCGKVFRG VAEFNEHRKS HVAAEPQPGP SRALENAAEK REQMEREAKP FECEECGKRF KKNAGLSQHL RVHSREKPF CEEGRSFKV NTHLFRHQKL

image



Background

Alternative Names Zinc finger protein 275

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZNF275 Protein
Catalog #: PEH6484

