


Product Name: Recombinant Human ZNF362 Protein
Catalog #: PEH6490



Summary

Name	ZN362_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5T0B9
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MPSQLDNLVL INKIKEQLMA EKIRPPHLPP TSASSQQPLL VPPAPAESSQ AVMSLPKLQQ VPGLHPQAVP QPDVALHARP ATSTVTGLGL STRTPSVSTS ESSAGAGTGT GTSTPSTPTT TSQSRLIASS PTLISGITSP PLLDSIKTIQ GHLLGPPKS ERGRKKIKAE NPGGPPVLVV PYPILASGET

image



Background

Alternative Names Zinc finger protein 362

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZNF362 Protein
Catalog #: PEH6490

