


Product Name: Recombinant Human ZNF615 Protein
Catalog #: PEH6491



Summary

Name	ZN615_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N8J6
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDVAVDFTWE EWQFLSPAQK DLYRDVMLEN YSNLVAVGYQ ASKPDALSKL ERGEETCTTE DEIYSRICSE IRKIDDPLQH HLQNQSIQKS VKQCHEQNMF GNIVNQKNGH FLLKQDCDTF DLHEKPLKSN LSFENQKRSS GLKNSAEFNR DGKSLFHANH KQFYTEMKFP AIAKPINKSQ FIKQQRTHNI ENAHVCSECG KAFLKLSQFI DHQRVHTGEK

image



Background

Alternative Names Zinc finger protein 615

Note

Product Name: Recombinant Human ZNF615 Protein
Catalog #: PEH6491



For Research Use Only , Not for Diagnostic Use