


Summary

Name	ZN618_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5T7W0
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RSQKSTKVEG PEPVPAEASL SAEQGTMTTEV KVKTELPDDY IQEVIWQGEA KEEKAVSKD GTSDVPAEIC VVIGGVRNQQ TLDGKAPEGS PHGGSVRSRY SGTWIFDQAL RYASGSYECG ICGKKYKYNN CFQTHVRAHR DTEATSGEGA SQSNNFRYTC DICGKKYKYY SCFQEHRDLH AVDVFSVEGA PENRADPFDQ GVVATDEVKE EPPEPFQKIG PKTGNYTCEF CGKQYKYYP YQEHVALHAP ISTAPGWEPP DDPDTGSECS HPEVSPSPRF VAAKTQTNQS GKKAPASVVR CATLLHRTTP ATQTQTFRTP NSGSPASKAT AAESAFSRRV EGKAQNHFEE

image



Background

Alternative Names Zinc finger protein 618

Product Name: Recombinant Human ZNF618 Protein
Catalog #: PEH6492



Note

For Research Use Only , Not for Diagnostic Use