


Product Name: Recombinant Human ZNF840P Protein
Catalog #: PEH6497



Summary

Name	ZN840_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A6NDX5
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RFNNNSNLNK HKKIHTGEKH FVCNQCGKAF SLNSKLSRHQ RTHNKKENSS KSVSNLNKHQ KTHAGEKPPF CNECKKAFAQ RMDLARHQQI

image



Background

Alternative Names Putative zinc finger protein 840 (Zinc finger protein 840 pseudogene)

Note

For Research Use Only , Not for Diagnostic Use